

697 | 801 L

■ ROUND long neck



Head Diameter 1/10mm	012	014	016	017	021	023	032	039
Head Length mm	1.2	1.4	1.6	1.7	2.1	2.3	3.2	3.9
ISO	order #							
806 315 697 504...	697-...XF				■021	■023		
806 315 697 514...	697-...F		■014		■021	■023		
806 315 697 524...	697-...M	■012	■014	■016	■017	■021	■023	■032 ■039
806 315 697 534...	697-...C	■012	■014	■016	■017	■021	■023	■032 ■039
806 315 697 544...	697-...XC				■017	■021	■023	

E164

Pointed cone w/rounded coated tip



Head Diameter 1/10mm	011	015	018
Head Length mm	10.0	10.0	10.0
ISO	order #		
806 315 164 534...	E164-...C	■011	■015 ■018

E165

Pointed cone uncoated tip



Head Diameter 1/10mm	014	014L
Head Length mm	9.0	9.0
ISO	order #	
806 315 165 534 ...	E165-...C	■014 ■014L

E213 "DA ACCESS"

SAFE END



Head Diameter 1/10mm	015
L1: Head Length mm	9.0
L2: OAL mm	26.5
ISO	order #
E806 316 213 524...	E213-...M ■015

152

ENDO ACCESS SAFE END



Head Diameter 1/10mm	016
L1: Head Length L1mm	9.0
L2: OAL mm	23.0
ISO	order #
500 315 210 295016	152-016 016

190 uncoated tip

FLAME OGIVAL SAFE END



Head Diameter 1/10mm	018
Head Length mm	9.5
ISO	order #
E806 314 190 534...	190-...C I018

216

TAPERED SAFE END



Head Diameter 1/10mm	010
Head Length mm	3.0
ISO	order #
806 314 216 504...	216-...XF I010
806 314 216 514...	216-...F I010

218

TAPERED SAFE END



Head Diameter 1/10mm	016
Head Length mm	6.0
ISO	order #
806 314 218 504...	218-...XF I016
806 314 218 514...	218-...F I016

219


TAPERED SAFE END



Head Diameter 1/10mm	012	016
Head Length mm	8.0	8.0
ISO	order #	
806 315 219 524...	219-...M I012 I016	


220 uncoated tip

TAPERED SAFE END



Head Diameter 1/10mm	014	016
L1: Head Length L1mm	10.0	10.0
ISO	order #	
806 315 220 524...	220-...M I014 I016	

Diamond Grit



GUIDE PINS

508C Cylinder

GUIDE PIN



Head Diameter 1/10mm	018	021
L1: Head Length L1mm	8.0	8.0
ISO	order #	
806 315 508C 514...	508C-...F	■ 018 ■ 021
806 315 508C 524...	508C-...M	■ 018 ■ 021



508T Taper

GUIDE PIN



Head Diameter 1/10mm	021	026
L1: Head Length L1mm	8.0	8.0
ISO	order #	
806 315 508T 514...	508T-...F	■ 021 ■ 026
806 315 508T 524...	508T-...M	■ 021 ■ 026



TC TUNGSTEN CARBIDE SURGICAL

151-Z



Head Diameter 1/10mm	016	
L1: Head Length L1mm	11.0	
L2: OAL mm	23.0	
ISO	order #	
500 315 199 336 016	151-Z	016



151M-Z



Head Diameter 1/10mm	016	
L1: Head Length L1mm	11.0	
L2: OAL mm	25.0	
ISO	order #	
500 316 199 336 016	151M-Z	016



151L-Z



Head Diameter 1/10mm	016	
L1: Head Length L1mm	11.0	
L2: OAL mm	28.0	
ISO	order #	
500 317 199 336 016	151L-Z	016



161

LINDEMANN straight



Head Diameter 1/10mm	016	
L1: Head Length L1mm	10.0	
L2: OAL mm	25.0	
order #	FG.TC000161	016



162

LINDEMANN cross cut



Head Diameter 1/10mm	016	
L1: Head Length L1mm	10.0	
L2: OAL mm	25.0	
order #	FG.TC000162	016

