



Polishine®



DCP

Daisy Contour Polishine®

New generation of open shape durable polishers with flexible lamellae

Daisy™ Contour Polishine
works flawlessly and extremely fast on both anterior and posterior restorations.

Daisy™ delivers a superb gloss and last much longer than any polisher.

The lamellae easily adapt to the surface structure, significantly increasing the contact area to the material.

This ensures a constant air supply and reduces frictional heat arising.

Composite

2-step polisher for high shine on all composite materials.

Polishers specifically adapted to each processing step.

Best results even on very hard composites.



Order #	DCP14SF	DCP14HG
Dimensions ø X L1	14 X 1.6	14 X 1.6
Color	Pink	Light Green
Grit	■ M	■ F
Recommended RPM	3,000 - 8,000	

Zirconia, Alumina

Especially suitable for full contour zirconia restorations.

High diamond concentration for more efficiency.

Makes glazing unnecessary.



Order #	DZP14SF	DZP14HG
Dimensions ø X L1	14 X 1.6	14 X 1.6
Color	Green	Orange
Grit	■ M	■ F
Recommended RPM	7,000 - 10,000	

Lithium Disilicate

Especially suitable for silicate porcelain.

Diamond grit for high efficiency.

Polishing result matches glazing.



Order #	DLP14SF	DLP14HG
Dimensions ø X L1	14 X 1.6	14 X 1.6
Color	Brown	Yellow
Grit	■ M	■ F
Recommended RPM	7,000 - 10,000	



90 | Polishing System

Polishine®

PLZ



Zirconia Polishing System

Including Autoclavable Bur Block

KIT PLZ - Zirconia Polishing System

2 steps Set

1st - Silk Finish : Dark Gray

2nd - Hi Gloss Polishing - Pink/Orange

Diamond grits combination for outstanding polishing results on Zirconia Oxide and hard Ceramics.

Pack of: 3+3

Autoclavable 135°C

RA SS shanks

1st step rec. RPM: 7,000-12,000. Max. 20,000

2nd step rec RPM: 7,000-10,000. Max. 20,000


Order No. kit PLZ

■ +FREE BUR BLOCK

1st - Silk Finish

2nd - Hi-Gloss Polishing

Zirconia



Item #	PLZ 101SF	PLZ 510SF	PLZ 675SF	PLZ 101HG	PLZ 510HG	PLZ 675HG
Dimensions Ø X L1	10 x 0.7	5 x 10	6.5 x 9.5	10 x 0.7	5 x 10	6.5 x 9.5
Color	Dark Grey	Dark Grey	Dark Grey	Pink/Orange	Pink/Orange	Pink/Orange
Grit	■ M	■ M	■ M	■ F	■ F	■ F
Recommended RPM	7,000 - 12,000	7,000 - 12,000	7,000 - 12,000	7,000-10,000	7,000-10,000	7,000-10,000



Polishine®

PLL

Designed for two steps polishing of Lithium Disilicate (e.g. IPS e.max®) restorations.

Lithium Disilicate System Including Autoclavable Bur Block

KIT PLL - Lithium Disilicate System

2 Steps Set

1st - Silk Finish: Blue

2nd - Hi Gloss Polishing: Yellow

Zirconium Oxide (ZrO₂) ceramic material for finishing and high gloss polish of Lithium Disilicate

Pack of: 3+3

Autoclavable 135°C

RA SS shanks

1st step rec. RPM: 7,000-12,000. Max. 20,000


2nd step rec. RPM: 7,000-10,000. Max. 20,000



Order No. kit PLL

+FREE BUR BLOCK

Lithium Disilicate Polishers



	1st - Silk Finish			2nd - Hi-Gloss Polishing		
	disc	cone	cup	disc	cone	cup
Item #	PLL 101SF	PLL 510SF	PLL 675SF	PLL 101HG	PLL 510HG	PLL 675HG
Dimensions Ø X L1	10 x 0.7	5 x 10	6.5 x 9.5	10 x 0.7	5 x 10	6.5 x 9.5
Color	Turquoise	Turquoise	Turquoise	Yellow	Yellow	Yellow
Grit	■ M	■ M	■ M	■ F	■ F	■ F
Recommended RPM	7,000 - 12,000	7,000 - 12,000	7,000 - 12,000	7,000-10,000	7,000-10,000	7,000-10,000

* IPS e.max is a registered trademark of IVOCCLAR VIVADENT, Inc.