



100 Carbide

ISO structure 500 314 168 006 010
 Carbide Shank ISO Grit Head
 type Shape Size ϕ 1/10mm

FINE TRIMMING & FINISHING 30 Blades

FG.TC

	Item	order #	Head ϕ 1/10 mm	Head Length L1 mm
	FG.TC EF9UF Safe End Esthetic Finisher	FG.TC0EF9UF	014	9.0
	FG.TC9103UF Interproximal	FG.TC9103UF	012	3.6
	FG.TC 9214UF Pointed taper	FG.TC9214UF	010	5.6
	FG.TC 9406UF Egg	FG.TC9406UF	018	3.4
	FG.TC 9903UF Needle	FG.TC9903UF	012	5.5

FG CARBIDE - SURGICAL

ROUND

FG.TC Surgical

	Item	ISO	Order #	Head ϕ 1/10 mm	Head Length L1 mm
	FG.TC 1/2SL	500 316 001 010 006	FG.TC1/2SL	006	0.6
	FG.TC 1SL	500 316 001 010 008	FG.TC001SL	008	0.8
	FG.TC 2SL	500 316 001 010 010	FG.TC002SL	010	1.0
	FG.TC 4SL	500 316 001 010 014	FG.TC004SL	014	1.4
	FG.TC 6SL	500 316 001 010 018	FG.TC006SL	018	1.8
	FG.TC 8SL	500 316 001 010 023	FG.TC008SL	023	2.3

STRAIGHT FLAT END FISSURE Cross Cut

FG.TC Surgical

	Item	ISO	Order #	Head ϕ 1/10 mm	Head Length L1 mm
	FG.TC 557SL	500 316 107 007 010	FG.TC0557SL	010	3.8
	FG.TC 558SL	501 316 107 007 012	FG.TC0558SL	012	3.8

TAPERED FISSURE Cross Cut


FG.TC Surgical	Item	ISO	Order #	Head Ø 1/10 mm	Head Length L1 mm
	FG.TC 701SL	500 316 168 007 012	FG.TC0701SL	012	3.8
	FG.TC 702SL	501 316 168 007 016	FG.TC0702SL	016	4.3

ZEKRYA

FG.TC Surgical	Item	ISO	Order #	Head Ø 1/10 mm	Head Length L1 mm
	151-Z (23mm OAL)	500 315 199 336 016	151-Z	016	11.0
	151M-Z (25mm OAL)	500 316 199 336 016	151M-Z	016	11.0
	151L-Z (28mm OAL)	500 317 199 336 016	151L-Z	016	11.0

New

LINDEMANN

FG.TC Surgical	Item	Order #	Head Ø 1/10 mm	Head Length L1 mm
	FG.TC 161	FG.TC000161	016	10.0
	FG.TC 162	FG.TC000162	016	10.0

FG CARBIDE - ENDO ACCESS

FG.TC Endo	Item	ISO	Order #	Head Ø 1/10 mm	Head Length L1 mm
	152-016	500 315 210 295 016	152-16	016	9.0