





Zirconia Ceramic Cutters

- Preparation of protheses.
- Preparation of Crowns, bridges and veneers.
- Aggressive yet precise and effective cutting.
- Effective substance removal.
- Finishing and polishing dentures.
- Removal of cement and luting residues after denture placement.

	Order #	Dimention	Shape
	HPZC 243CC XC	14.5 X 6.0	Bullet
	HPZC 243SC C	14.5 X 6.0	Bullet
	HPZC 243SC F	14.5 X 6.0	Bullet
	Order #	Dimention	Shape
	HPZC 201CC XXC	14.5 X 6.0	Conical
	HPZC 201CC XF	14.5 X 6.0	Conical
	HPZC 201SC XC	13.0 X 4.5	Conical
	HPZC 175SC XC	17.0 X 2.3	Conical